


Decalcification Products



Super Decalcifier I - *Delicate*

Recommended for use with delicate specimens that require fast diagnostic work such as IHC, cytochemistry and special stain procedures.


- Save Time - Decalcify in as little as 3 hours!
- Ready-to-Use - for human and animal tissue
- Recommended for IHC, special stains and routine stains
- Enhanced nuclear detail
- Environmentally safe and biodegradable



Super Decalcifier II - *Heavy Duty*

Extremely effective and versatile decalcifier that can be used in specific lab routines. As with most acids, nuclei acids in the cell can become subject to ribonuclease digestion, resulting in a loss of basophilic properties.

- Save Time - Decalcify in as little as 3 hours!
- Ready-to-Use - decalcifier for human and animal tissues
- Enhanced nuclear detail
- Environmentally safe and biodegradable



Deli-Cal Block Solution

Our microtome soak, Deli-Cal Block Solution saves time and money by allowing to salvage problem blocks right at cutting stations. The quality and thoroughness of Deli-Cal Block Solution is dependent upon the original decalcification procedure and condition of the specimen.

- Saves knives and saves money
- Useful for Bone Marrows, Special Stains & IHC
- Enhances ribbon cutting
- Salvages undecalcified and under processed bone marrow biopsies at the cutting station
- Fast soaking procedure
- Environmentally safe and biodegradable



De-Calcify Block Solution

Original decalcification thoroughness and condition of the specimen is dependent on the original decalcification solution. Use De-Calcify Block Solution as a fast microtome soak that saves you time and money by allowing you to salvage problem blocks right at your cutting station.

- Saves knives and saves money
- Enhances ribbon cutting
- Useful for H&E and routine staining
- Fast soaking procedure
- Ready-to-use formulation
- Environmentally safe and biodegradable



Cat. # Sizes

Super Decalcification I - *Delicate*

24888 500ml, 1L, 1 case (6x1 liter)

Super Decalcification II - *Heavy Duty*

24887 500ml, 1L, 1 case (6x1 liter)

Deli-Cal Block Solution

24900 250ml, 500ml

De-Calcify Block Solution

24903 250ml, 500ml