

PolyTransport Buffer



Preserve Your Specimens “To Go”

Preserve specimens from surgery to the point of grossing for fixation prior to further processing.

Specially formulated transport medium designed to preserve tissue in a viable state during transport by maintaining osmolarity and pH for 24 to 48 hours.



Benefits and Advantages

- Preserves tissue integrity during transport
- Prevents proteins from denaturing or crosslinking
- Prevents bacterial growth
- Contains a cryoprotectant to prevent ice crystal formation during rapid freezing techniques
- Contains a membrane stabilizing agent to prevent lysis
- No washing of specimens necessary prior to placing in fixative or rapidly freezing

Size: 250ml, 500ml

Cat. #: 24311